

# PROBES

## Gas Probes

Due to a broad range of gas probes with standardized connection forms (DIN, BS, NFS, AGA), Schön connection tubes can be used worldwide for direct extraction from the central gas supply system. In addition to high safety technology, the gas probes distinguish themselves by comfortable handling.

### Medical Gas Probes Function

Medical gas probes are designed to be used on all brands of gas outlets. The most common types are shown below. The advantage of the medical gas probes is the snap-in/snap-out connection it offers.

### Features

- Fast and easy to connect and disconnect
- Labeled and indexed for specific gas services
- Available in NPT and hose styles



BS 5682:2015 Standard Probe



DIN 13260 Standard Probe



NF S 90-116 Standard Probe



AGA Standard Probe

## DIN , BS, NF-S, AGA Standard Medical Gas Probe



DIN 13260-2



BS 5682:2015



AGA

NFS 90-116

